Film-Coated Tablet



Diosmin + Hesperidin

### Description

Pilestop® is a fixed dose combination of Diosmin and Hesperidin. It is a phlebotonic drug and a vascular protecting agent. It increases venous tone by prolonging the vasoconstrictor effect of noradrenaline on the vessel wall, and reduces blood venous stasis. It reduces the inflammatory process at the microcirculatory level, by inhibiting the adhesion, migration, and activation of leukocytes. This reduces capillary hyperpermeability and increases capillary resistance. It increases lymphatic contractions and improves lymphatic drainage.

### Composition

Each film-coated tablet contains micronised Diosmin BP 450 mg and Hesperidin INN 50 mg.

#### Indications

- · Acute hemorrhoidal attacks.
- Chronic hemorrhoidal disease.
- Chronic venous insufficiency.
- Lymphedema

### **Dosage & Administration**

Tablet should be taken at meal times.

Disease condition	Morning	Noon	Evening	Duration of treatment
Acute hemorrhoidal attack	2 tablets	2 tablets	2 tablets	First four days
	2 tablets	-	2 tablets	Next Three days
Maintenance dose (If required)	1 tablet	-	1 tablet	Approx. two months
Chronic hemorrhoidal disease	1 tablet	-	1 tablet	Approx. two to three months
Chronic venous insufficiency (CVI)	1 tablet	-	1 tablet	Two months
Lymphedema	1 tablet	-	1 tablet	Six months

### OR AS DIRECTED BY THE PHYSICIAN.

#### Side Effects

Some cases of minor gastrointestinal and autonomic disorders have been reported, but these never required cessation of treatment.

# Contraindications

Data is not available yet.

## Use in Pregnancy and Lactation

Pregnancy: Experimental studies in animal have not demonstrated any teratogenic effects and no harmful effects have been reported in man to date.

Lactation: In the absence of data concerning the diffusion into breast milk, breast-feeding is not recommended during treatment.

## Precaution

Data is not available yet.

# **Drug Interaction**

Data is not available yet.

### Supply

Each box containing 2 x 10's tablets in blister pack.

Store in a cool & dry place protected from light. Keep all medicines out of reach of children.

\* Further information is available on request.

Manufactured by:

The ACME Laboratories Ltd.

Dhamrai, Dhaka, Bangladesh

07 1158/0

